

## Product Information

FreezeMe One, Cryopreservation Medium with Fetal Bovine Serum

Cat. No.: FM1-F (50 ml)

---

### Product Description

FreezeMe One is a cryoprotective medium that contains 10 % DMSO and 20 % FBS. The serum used is extensively tested to protect cells during cell preservation.

### Applications

- Cryopreservation of a wide range of cell types with high viability
- Ready-to-use solution

### Product Specifications

Sterility	Tested
Storage	≤-15°C. FreezeMe One is a light sensitive solution. It should be protected from light during shipping and storage.

---

### Freezing Protocol

---

Before cryopreservation cells should be checked for contamination. FreezeMe One can be used with any standard freezing protocol.

#### Cryopreservation of Suspension Cultures

- Count the number of viable cells to be cryopreserved. Cells should be in mid-log phase of growth. Centrifuge for 5 min to pellet cells (200 to 400 g). Remove the supernatant down to the smallest volume without disturbing the cells.
- Resuspend cells in pre-cooled (+4°C to +8°C) FreezeMe One to a concentration of  $5 \times 10^6$  to  $10^7$  cells/ml.
- Aliquot into cryogenic storage vials. Place vials at +4°C and start the freezing procedure within 5 min. Cells are frozen slowly at +1°C/min (by programmable coolers or by placing vials in an insulated box in a -70°C to -90°C freezer).
- Then transfer storage vials to liquid nitrogen storage.

#### Cryopreservation of Adherent Cultures

- Detach cells from the substrate with a gentle dissociating agent. Especially with sensitive cells use Accutase (Cat. No. ACC-1B) to avoid cell damage. Inactivate dissociating agent if necessary.
  - Resuspend the detached cells in complete growth medium and establish the viable cell count.
  - Centrifuge for 5 min to pellet cells (200 to 400 g). Remove the supernatant down to the smallest volume without disturbing the cells.
  - Resuspend cells in pre-cooled (+4°C to +8°C) FreezeMe One to a concentration of  $5 \times 10^6$  to  $10^7$  cells/ml.
  - Aliquot into cryogenic storage vials. Place vials at +4°C and start the freezing procedure within 5 min. Cells are frozen slowly at +1°C/min (by programmable coolers or by placing vials in an insulated box in a -70°C to -90°C freezer).
  - Then transfer storage vials to liquid nitrogen storage.
-

## Product Information

FreezeMe One, Cryopreservation Medium with Fetal Bovine Serum

Cat. No.: FM1-F (50 ml)

---

### Thawing of Cryopreserved Cells

---

Cryopreserved cells can be thawed by the following procedures:

#### Centrifugation

- Remove cells from storage and thaw quickly in a +37°C water bath. Capricorn Scientific recommends eye protection by using approved safety goggles. We also suggest the use of safety gloves to protect uncovered skin.
- Place 1 to 2 ml of thawed cells in ~25 ml of complete growth medium. Mix cell suspension gently.
- Centrifuge the cells at ~80 g for 2 to 3 min.
- Check clarity of the supernatant and visibility of a consolidated cell pellet. Discard supernatant without disturbing the cells.
- Gently resuspend the cells in complete growth medium and perform a viable cell count.
- Plate the cells. Cell inoculum should be at least  $3 \times 10^5$  viable cells/ml.

#### Direct plating

- Remove cells from storage and thaw quickly in a +37°C water bath. Capricorn Scientific recommends eye protection by using approved safety goggles. We also suggest the use of safety gloves to protect uncovered skin.
- Plate cells directly with complete growth medium. Use 10 to 20 ml of complete medium per 1 ml of frozen cells. Cell inoculum should be at least  $3 \times 10^5$  cells/ml.
- Culture cells for 12 to 24 h. Replace medium with fresh complete growth medium to remove cryopreservative.

We recommend thawing procedure 1, especially when handling sensitive cells.

#### Precautions and Disclaimer

This product is for research use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

**Caution:** FreezeMe One contains Dimethyl Sulfoxide (DMSO). Do not breathe gas/fumes/vapour/spray. Avoid contact with eyes and skin. Irritant to eyes, respiratory system and skin. S23 S24/25.

#### Help needed?

If you have any further questions regarding this product please do not hesitate to contact our cell culture experts:

Email: [techservice@capricorn-scientific.com](mailto:techservice@capricorn-scientific.com)

Phone: +49 6424 9446427

Accutase is a registered trademark of Innovative Cell Technologies, Inc.