

ELECTROPHORESIS

PROTEAN® II xi and XL Cells

- Large format vertical electrophoresis systems
- PROTEAN II xi format for 16 x 16 cm or 16 x 20 cm gels
- PROTEAN II XL format for 2-D gels with 17 cm IPG strips
- Cells accommodating 2 to 4 gels; multi-cells accommodating up to 6 gels
- Enhanced resolution over smaller formats

PROTEAN II xi and XL Cells and Multi-Cells for Large Format Electrophoresis

Large Format Vertical Electrophoresis

The PROTEAN II xi and PROTEAN II XL cells and multi-cells are large format vertical electrophoresis cells that offer the versatility to perform a wide variety of separations, including SDS- and native PAGE, 1-D and 2-D PAGE, preparative and gradient PAGE, and even high-resolution nucleic acid separations on agarose gels. For 1-D electrophoresis, use the PROTEAN II xi cell or xi multi-cell to run 16 x 16 cm (precast) or 16 x 20 cm (handcast) gels. For 2-D analysis, use the PROTEAN II XL cell or XL multi-cell to run 18.3 x 19.3 cm (precast) or 18.5 x 20 cm (handcast) second-dimension gels that accommodate 17 cm immobilized pH gradient (IPG) strips.

Whatever the application, large gel formats for vertical electrophoresis offer increased separation area, enable analysis of larger sample loads, and ultimately provide superior resolution over smaller gel formats.

Easy to Assemble

- Single screw clamps minimize assembly time
- Casting seal eliminates messy grease and agarose leaks when hand casting gels

Versatile

- Run 2 to 4 slab gels simultaneously with the PROTEAN II xi and XL cells or up to 6 gels with the multi-cells
- For second-dimension runs using 17 cm IPG strips, expand the PROTEAN II xi cell to the PROTEAN II XL format with the PROTEAN II IPG conversion kit



PROTEAN II xi 2-D Precast Gel System

- Use with handcast or PROTEAN II xi or XL Ready Gel® precast gels, both single percentage or gradient. All are available with IPG well for 2-D electrophoresis and in 15- or 20-well formats for 1-D electrophoresis. Gels are available for both the PROTEAN II xi and PROTEAN II XL formats
- Utilize accessories available for hand casting gels
- The PROTEAN II xi 2-D cell is available for running tube gels in the first dimension and slab gels in the second dimension

Reliable

- Central cooling core maximizes resolving power and gel-to-gel reproducibility
- Buffer recirculation maintains uniform pH for better separation and reproducibility
- Raised upper electrode prevents arcing
- Safety lid prevents accidental reversal of polarity
- High-quality materials offer durability

BIO-RAD

Ordering Information

Catalog #	Description
PROTEAN II xi Cells* for 16 x 16 cm and 16 x 20 cm Gel Sizes	
165-1801	PROTEAN II xi Cell , 16 cm, without spacers and combs
165-1802	PROTEAN II xi Cell , 16 cm, 1.5 mm spacers (4), 15-well combs (2)
165-1803	PROTEAN II xi Cell , 16 cm, 1.0 mm spacers (4), 15-well combs (2)
165-1804	PROTEAN II xi Cell , 16 cm, 0.75 mm spacers (4), 15-well combs (2)
165-1811	PROTEAN II xi Cell , 20 cm, without spacers and combs
165-1812	PROTEAN II xi Cell , 20 cm, 1.5 mm spacers (4), 15-well combs (2)
165-1813	PROTEAN II xi Cell , 20 cm, 1.0 mm spacers (4), 15-well combs (2)
165-1814	PROTEAN II xi Cell , 20 cm, 0.75 mm spacers (4), 15-well combs (2)
PROTEAN II XL Cells, Wide Format, Compatible With ReadyStrip™ IPG Strips, for 18.3 x 20 cm Gel Size	
165-3188	PROTEAN II XL Cell, wide format , 1.0 mm, includes PROTEAN II xi basic unit (165-1834) and 1.0 mm IPG conversion kit (165-3183)
165-3189	PROTEAN II XL Cell, wide format , 1.5 mm, includes PROTEAN II xi basic unit (165-1834) and 1.5 mm IPG conversion kit (165-3186)
165-3190	PROTEAN II XL Cell, wide format , 2.0 mm, includes PROTEAN II xi basic unit (165-1834) and 2.0 mm IPG conversion kit (165-3184)
PROTEAN II IPG Conversion Kits for 2-D Electrophoresis (to Convert xi to XL)	
165-3183	PROTEAN II xi Cell IPG Conversion Kit , 1.0 mm, 18.5 x 20 cm, includes 2 sets of IPG clamps, 2 sets of 20 x 20 cm glass plates, 4 IPG spacers, 2 IPG 2-D combs, 2 IPG central cooling core gaskets, 2 casting stand gaskets, alignment card
165-3186	PROTEAN II xi Cell IPG Conversion Kit , 1.5 mm, 18.5 x 20 cm
165-3184	PROTEAN II xi Cell IPG Conversion Kit , 2.0 mm, 18.5 x 20 cm
165-1834	PROTEAN II xi Basic Unit With Casting Stand , includes central cooling core, lower buffer chamber, lid with cables, leveling bubble; combine with an IPG conversion kit for a complete 18.5 cm wide format system
PROTEAN II xi Multi-Cell	
165-1951	PROTEAN II xi Multi-Cell , includes 3 central cooling cores with gaskets, buffer tank, lid with power cables, 1 upper buffer dam, PROTEAN II xi multi-casting chamber with accessories, leveling bubble, instructions
165-1956	PROTEAN II xi Multi-Cell 2-D Conversion Kit , for proper cooling, includes 2 cooling coils and manifold (required for 2-D electrophoresis applications)
PROTEAN II XL Multi-Cell, Wide Format, Compatible With ReadyStrip IPG Strips	
165-3176	PROTEAN II XL Multi-Cell , wide format, 1.0 mm
165-3177	PROTEAN II XL Multi-Cell , wide format, 1.5 mm
165-3178	PROTEAN II XL Multi-Cell , wide format, 2.0 mm

Catalog #	Description
PROTEAN II xi 2-D Cells**	
165-1931	PROTEAN II xi 2-D Cell , 1.0 mm, 16 cm
165-1932	PROTEAN II xi 2-D Cell , 1.5 mm, 16 cm
165-1933	PROTEAN II xi 2-D Cell , 1.0 mm, 20 cm
165-1934	PROTEAN II xi 2-D Cell , 1.5 mm, 20 cm
PROTEAN II XL Ready Gel Precast Gels, IPG Well	
161-1450	10% Tris-HCl Gel , 1.0 mm thick, 18.3 cm x 19.3 cm, IPG well
161-1451	12% Tris-HCl Gel , 1.0 mm thick, 18.3 cm x 19.3 cm, IPG well
161-1452	10–20% Tris-HCl Gel , 1.0 mm thick, 18.3 cm x 19.3 cm, IPG well
161-1453	8–16% Tris-HCl Gel , 1.0 mm thick, 18.3 cm x 19.3 cm, IPG well
PROTEAN II xi Ready Gel Precast Gels, 15-Well	
161-1454	10% Tris-HCl Gel , with 4% stacking gel, 1.0 mm thick, 16 x 16 cm, 15-well
161-1455	12% Tris-HCl Gel , with 4% stacking gel, 1.0 mm thick, 16 x 16 cm, 15-well
161-1456	10–20% Tris-HCl Gel , with 4% stacking gel, 1.0 mm thick, 16 x 16 cm, 15-well
161-1457	8–16% Tris-HCl Gel , with 4% stacking gel, 1.0 mm thick, 16 x 16 cm, 15-well
PROTEAN II xi Ready Gel Precast Gels, 20-Well	
161-1458	10% Tris-HCl Gel , with 4% stacking gel, 1.5 mm thick, 16 x 16 cm, 20-well
161-1459	12% Tris-HCl Gel , with 4% stacking gel, 1.5 mm thick, 16 x 16 cm, 20-well
161-1460	10–20% Tris-HCl Gel , with 4% stacking gel, 1.5 mm thick, 16 x 16 cm, 20-well
161-1461	8–16% Tris-HCl Gel , with 4% stacking gel, 1.5 mm thick, 16 x 16 cm, 20-well

* Each PROTEAN II xi cell includes a central cooling core with gaskets, lower buffer chamber, lid with cables, 2 sets of glass plates, 4 sandwich clamps, upper buffer dam, casting stand with gaskets, leveling bubble, and instructions. The 16 cm cells include clamps for running 16 cm gels; the 20 cm cells include clamps for 20 cm gels. Cells without spacers and combs (catalog #165-1801 and #165-1811) allow you to select the spacer thickness and comb that will best fit your application.

** Each PROTEAN II xi 2-D cell includes a central cooling core with gaskets, lower buffer chamber, lid with cables, 2 sets of glass plates (with beveled inner plates), 4 sandwich clamps, twenty-four 180 mm long glass tubes (tube diameter = spacer thickness), 2 tube gel adaptors, 16 grommets, 16 stoppers, two 2-D combs, 4 spacers, upper buffer dam, casting stand with gaskets, leveling bubble, and instructions. Sandwich clamps are sized to fit the gel length appropriate for the cell (16 cm or 20 cm). 1.0 mm and 1.5 mm indicate the thickness of spacers and combs included with the cell. Run both tube gels for first dimension electrophoresis and slab gels for second dimension electrophoresis in the same cell.



**Bio-Rad
Laboratories, Inc.**

Life Science
Group

Web site www.bio-rad.com USA 800 4BIORAD Australia 61 02 9914 2800 Austria 01 877 89 01 Belgium 09 385 55 11 Brazil 55 21 3237 9400
Canada 905 364 3435 China 86 21 6426 0808 Czech Republic 420 241 430 532 Denmark 44 52 10 00 Finland 09 804 22 00 France 01 47 95 69 65
Germany 089 318 84 0 Greece 30 210 777 4396 Hong Kong 852 2789 3300 Hungary 36 1 455 8800 India 91 124 4029300 Israel 03 963 6050
Italy 39 02 216091 Japan 03 6361 7000 Korea 82 2 3473 4460 Mexico 52 555 488 7670 The Netherlands 0318 540666 New Zealand 0508 805 500
Norway 23 38 41 30 Poland 48 22 331 99 99 Portugal 351 21 472 7700 Russia 7 495 721 14 04 Singapore 65 6415 3188 South Africa 27 861 246 723
Spain 34 91 590 5200 Sweden 08 555 12700 Switzerland 061 717 95 55 Taiwan 886 2 2578 7189 United Kingdom 020 8328 2000